



# *UV Activator*



Attractive to **water**,  
Attractive to **patients!**

A 20-second Hydrophilic Surface  
Treatment Solution!

**DIO** IMPLANT




## Intuitive User Interface

### Power

White – Ready to activate  
Red – Stand by mode

### UV Indicator

Light on – UV cycle active  
Light off – UV cycle complete



### Door

Light on - Door open or ajar  
Light off - Door closed

### Dimensions

13.8"(h) X 10.3"(w) X 7.1"(d)  
Power consumption : 150W maximum  
Weight : 6kg  
Recommended UV exposure : 20 seconds  
Activates one fixture per cycle

### \* Two color options



Rose Pink

Sky Blue

Activate in just 20 seconds with full surround UV light

## 20 miraculous seconds

Experience maximal hydrophilic qualities!

	Wet	UV-A Wavelength	Dry DIO UV Activator
Activation method	Activated during the manufacturing process, delivered in solution-filled packaging	UV activated before use	UV activated before use
Activation time		20 hours ↑	Only 20sec

← 3,600 times faster! →

## Is UV Active right for you?

Consider these:

For elderly patients with weak or insufficient bone

**INDICATED!**

For patients with a busy schedule

**INDICATED!**

For immediate loading or large post-extraction socket bony defect cases

**INDICATED!**

# Up your game with UV!

- A truly hydrophilic surface with only 20 seconds of UV Activation
- Bathed in full surround UV light ensuring maximum activation
- Intuitive user interface, alluring design



# ***UV Activator***

63mm x 17mm

\*본 기기는 의료기기가 아님  
UV Activator는 임플란트와 같은 의료기기 등을 자외선 파장 처리하여  
친수성을 도모하는 기기입니다.  
(방송통신기자재등의 적합 인증 : R-C-DIO-UVACT02)